



Product Overview

Track End Caps are designed to provide a secure, professional finish for patient-lifting track installations. Installed at both ends of the track, they close the track profile, protect internal components, and help maintain the integrity and appearance of the completed system.



Track End Cap Dimensions (mm/inches)			
Dimensions (Ref to drawing)	Single Track End Cap	Mid-Duty Track End Cap	Heavy-Duty Track End Cap
A	2.68 in	5.82 in	7.08 in
B	2.28 in	2.28 in	2.28 in
C	0.78 in	0.78 in	0.78 in
D	0.92 in	0.92 in	0.92 in
E	1.38 in	1.38 in	1.38 in
F	0.4 in	0.4 in	0.4 in
G	0.24 in	0.24 in	0.24 in



- Two End Caps are required at both ends of the system (one at each end).
- The End Caps are used to close the installation, creating a complete system.
- End caps are used in combination with the safety bolt to secure them in place.

Please refer to the installation manual for the correct installation method.